
(DE) Digital Thermometer #2136 / #2137

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(EN) Digital Thermometer #2136 / #2137

Note: The exterior of each model has a little difference.

Please read the instructions carefully before using the thermometer for the first time, and keep these in a safe place.

Correct use
The fever thermometer is designed exclusively for use in the veterinary sector and not for use on humans.

Use
Before using, please disinfect the probe at first. To switch on, press the ON/OFF button next to the display; a short beep will sound, indicating that the thermometer is operational. When the letters „Lo" and a flashing „°C" display, the thermometer is now ready for use. If the ambient temperature is below 32 °C, then „Lo°C" will appear on the LCD and if it is more than 42°C, then „HI°C" will appear on the LCD. As soon as the constant temperature value is reached, a beep will sound four times, and the „°C" symbol will stop flashing. The highest temperature measured appears on the LCD. When the measurement is completed, please switch the thermometer off by pressing the ON/OFF button. After the temperature has been displayed, the thermometer will shut off automatically after 10 minutes.

Memory function
Switch the thermometer off, then press the ON/OFF button for about 2 seconds. The last measured value with „°C" will appear automatically on the LCD.

Methods of measuring temperature
In the rectum (rectal)
This is the most accurate method from a medical point of view. The thermometer tip is inserted carfully into the rectum for a maximum of 2 cm. The usual measuring time is approximately 9 to 11 seconds.

In the mouth (oral)
There are different heat zones in the mouth. As a general rule, the oral temperature is 0.3 °C to 0.8 °C lower than the rectal temperature. The usual measuring time is approximately 15 to 19 seconds.

Note: We strongly recommend the rectal method as the most accurate method for identifying the basal temperature and advise you to extend the measuring time by 3 minutes after the beep.

Cleaning and disinfection
The best way to clean the thermometer tip is by applying a disinfectant (e.g. 70 % ethyl alcohol) with a damp cloth.

Safety precautions
- Do not allow the device to come into contact with hot water.
- Do not expose to high temperateres or direct sunlight.
- Do not drop the thermometer. It is neither shock-proof nor impact-resistant.
- Do not bend or open the device (except the battery compartment).
- Do not clean with thinners, petrol or benzenel. Only clean with water or disinfectant.
- Do not immerse the Waterproof thermometers in liquid for longer than 30 minutes.
- The thermometer contains small parts (battery, battery compartment) which can be swallowed by children. For this reason, do not leave the thermometer unattended in the hands of children.
- Avoid bending the Non-Flexible thermometer tip. For the Flexible thermometer, avoid bending the tip by more than 45 degrees.
- If the ambient temperature is over 35 °C, dip the thermometer tip in cold water for approx. 5 to 10 seconds prior to measuring the temperature.
- A persistent fever requires medical attention - please consult a veterinarian!

Battery replacement
The battery is empty and needs replacing when the "▬" or "▭" battery symbol appears on the right of the LCD. Remove the battery cover and replace it with a battery (preferably non-mercury) of the same type. Please note: The "+ " sign up and "-" sign down.

Technical data

Type:	maximum thermometer
Measurement range:	(32.0 ~ 42.0) °C
Measurement accuacy:	+/- 0.1 °C (35.5 °C ~ 42.0 °C) +/- 0.2 °C (for the remaining measurement range)
Storage/transportation temperature:	(-25 ~ 55) °C, ⩽ 95 % RH
Ambient temperature during use:	(5 ~ 35) °C ⩽ 80 % RH
Min Scale:	0.1 °C
Battery type:	Alkallmanganese battery, type LR41, 1.5 V, service life minimum 100 hours under continuous operation.
Weight:	Approx. 13 g

Calibration check
This thermometer is initially calibrated at the time of manufacture. If this thermometer is used according to the operation instruction, periodic re-adjustment is not required.

CE-/UKCA-conformity declaration
Albert KERBL GmbH hereby declares that the product / device described in these instructions complies with the fundamental requirements and other relevant stipulations and regulations. The CE-/UCKA-mark confirms compliance with the Directives of the European Union or the relevant UK legislation.

Electronic scrap
Disposing of this device after its service life is the responsibility of the operator. Please consult the valid national regulations. The device must not be disposed of in household waste. In accordance with the stipulations of the EU Directive on the Disposal of Electrical and Electronic Devices, the device can be disposed of free of charge at the local waste collection or recycling centre. Alternatively, it can be returned to retailers who offer a collection service. The proper disposal helps to ensure environmental protection and prevents any adverse effects on human health and the environment.

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Albert Kerbl GmbH Felizenzell 9
Kerbl UK Ltd Lands End Way